

Conjunto PH-HPS



El conjunto Photonic HPS está compuesto por un foco LED de alta potencia con dos brazos y por los adaptadores MS35/36B(E)

Este foco se utiliza para unirse directamente a los soportes del microscopio, como los soportes Dino-Lite MS35/36B(E). Gracias a que se monta sobre el soporte, es perfecto para espacios pequeños. La luminosidad se puede controlar con la unidad de control compacta. Además, puede encender y apagar los dos brazos de manera individual.

Los brazos flexibles pueden adaptarse fácilmente para proporcionar el mejor ángulo de iluminación a su muestra.

* El resto de los artículos de Photonic está disponible en toda Dino-Lite Europe.

Información adicional

- Accessories range: Light and control
- Dimensions: 55cm long, each arm
- Weight: 650 g
- ESD Safe: Yes
- Compatible with: All Dino-Lite models, except high magnification models
- Package contents: Hi-Power Spot 2-arms, Controller unit, power supply with EU & UK plug, MS35/36B(E) adapter | ESD Safe add-ons not included
- Price range: €450,00 - €550,00

- **Details**

- Accessory range: Light and control
- Price range: €450,00 - €550,00
- ESD-Safe: Yes
- Dimensions: 55cm long, each arm
- Weight: 650 g
- Compatible with: All Dino-Lite models, except high magnification models
- Package contents: Hi-Power Spot 2-arms, Controller unit, power supply with EU & UK plug, MS35/36B(E) adapter | ESD Safe add-ons not included
- Warranty: 2 year European warranty
- Accessory type: Light and control
- Price range (excl. VAT): €450,00 - €550,00
- ESD-safe: Yes

PDF generated on: 24-10-2018