

Interruptor de pie SW-F1



El SW-F1 le permite tomar una fotografía con una simple presión del pie utilizando soluciones de microscopía digital con los productos Dino-Lite o Dino-Eye. El "interruptor de pie" sólo es compatible con DinoCapture 2.0 para PC y DinoXcope para los ordenadores Macintosh. No se necesita software adicional, basta con conectar el dispositivo al puerto USB y está listo para usar.

Información adicional

- Accessories range: Light and control
- Dimensions: 10cm (W) x 7cm (L) x 3cm (H)
- Weight: 300 g
- ESD Safe: Yes
- Compatible with: All USB models
- Package contents: 1 foot pedal
- Price range: €0 - €50,00

• Details

- Accessory range: Light and control
- Price range: €0 - €50,00
- ESD-Safe: Yes
- Dimensions: 10cm (W) x 7cm (L) x 3cm (H)
- Weight: 300 g
- Compatible with: All USB models
- Package contents: 1 foot pedal

- Warranty: 1 year European warranty
- Accessory type: Light and control
- Price range (excl. VAT): €0 - €50,00
- ESD-safe: Yes

PDF generated on: 24-10-2018