

KM-01



El KM-01 se puede usar para controlar la rueda de enfoque de Dino-Lite sin tocar el dispositivo. Esto resulta especialmente útil cuando hay que aislar el Dino-Lite mientras se utiliza.

Los controles táctiles situados en la robusta caja de control del KM-01 proporcionan una respuesta inmediata, pero el KM-01 también se puede controlar con el software DinoCapture 2.0. El KM-01 se puede montar sobre los soportes MS35B/MS36B/RK-10A de Dino-Lite, o sobre soportes de terceras partes con un mástil de 10 mm de diámetro.

El KM-01 es compatible con determinados modelos de Dino-Lite. Rogamos consulte a su distribuidor local o a Dino-Lite Europe para obtener información adicional.

Información adicional

- Accessories range: Light and control
- Dimensions: 77.14mm (L) 65.55mm (W) 98mm (H)
- Weight: 340 g (Motor) 110 g (Control Box)
- ESD Safe: No
- Compatible with: Please inquire with your local Dino-Lite sales representative for compatibility information
- Package contents: Motor, Control Box, C-Clip, Set Screws, AC Adapter, USB Cable, Male/Female 2 Pin Circular DIN Connector, Software CD
- Price range: €250,00 - €350,00

- **Details**

- Accessory range: Light and control
- Price range: €250,00 - €350,00
- ESD-Safe: No
- Dimensions: 77.14mm (L) 65.55mm (W) 98mm (H)
- Weight: 340 g (Motor) 110 g (Control Box)
- Compatible with: Please inquire with your local Dino-Lite sales representative for compatibility information
- Package contents: Motor, Control Box, C-Clip, Set Screws, AC Adapter, USB Cable, Male/Female 2 Pin Circular DIN Connector, Software CD
- Warranty: 1 year European warranty
- Accessory type: Light and control
- Price range (excl. VAT): €250,00 - €350,00
- ESD-safe: No

PDF generated on: 24-10-2018