

MS33W



Base metálica redonda con brazo flexible y soporte para el Dino-Lite. El MS33W puede modificar la altura (aprox. 15 mm) usando un mando giratorio en la base.

El soporte de escritorio MS33W sostiene el Dino-Lite con un brazo flexible con cuello de cisne que se puede doblar y girar en muchas posiciones. Una funcionalidad adicional del MS33W es el mando de ajuste preciso situado en la base del soporte. Utilizando este mando, usted puede enfocar el Dino-Lite en situaciones en las que puede resultar muy difícil mover el disco de ajuste en el propio microscopio.

Información adicional

- Accessories range: Basic stands
- Dimensions: 13.3cm (W) x 10.1cm (L) x 30.5cm (H)
- Weight: 500 g
- ESD Safe: No
- Compatible with: All Dino-Lite models
- Package contents: Stand, HL2 holder
- Price range: €0 - €50,00

• Details

- Accessory range: Basic stands
- Price range: €0 - €50,00

- ESD-Safe: No
- Dimensions: 13.3cm (W) x 10.1cm (L) x 30.5cm (H)
- Weight: 500 g
- Compatible with: All Dino-Lite models
- Package contents: Stand, HL2 holder
- Warranty: 1 year European warranty
- Accessory type: Basic stands
- Price range (excl. VAT): €0 - €50,00
- ESD-safe: No

PDF generated on: 24-10-2018