

MS52BA2



Base de mesa de alto rendimiento para el MS52SB. El MS52BA2 es una combinación del brazo flexible articulado MS52B y un soporte base de escritorio que sostiene al MS52B encima de una mesa.

El soporte base de escritorio está diseñado para situaciones en las que no se pueda utilizar la abrazadera estándar para sujetar el MS52B en una superficie estable. El brazo flexible se mueve y ajusta apretando los mandos y las abrazaderas para posicionar el brazo en el área que usted desea examinar. Dos de los ejes del brazo flexible se pueden girar 360 grados y se pueden ampliar para obtener un mayor alcance. El brazo flexible MS52BA2 está fabricado con aluminio reciclado duradero y ligero y posee la capacidad de ser plegado, lo que hace que el soporte MS52BA2 sea portátil y adecuado para viajes.

Información adicional

- Accessories range: Professional stands
- Dimensions: 35.5cm (W) x 16.5cm (L) x 5cm (H)
- Weight: 4 kg
- ESD Safe: No
- Compatible with: All Dino-Lite models
- Package contents: Stand, table base, HD-P1 holder
- Price range: €250,00 - €350,00

- **Details**

- Accessory range: Professional stands
- Price range: €250,00 - €350,00
- ESD-Safe: No
- Dimensions: 35.5cm (W) x 16.5cm (L) x 5cm (H)
- Weight: 4 kg
- Compatible with: All Dino-Lite models
- Package contents: Stand, table base, HD-P1 holder
- Warranty: 1 year European warranty
- Accessory type: Professional stands
- Price range (excl. VAT): €250,00 - €350,00
- ESD-safe: No

PDF generated on: 24-10-2018