

EarScope Basic



El EarScope Basic (MEDL3E) es una solución compacta y segura para realizar investigación de oídos. La cápsula es desmontable para fines higiénicos. El MEDL3E tiene una resolución de 640x480 y un aumento de 20 ~ 30x. Este modelo se ofrece con espéculos desechables de 3, 4 y 5 mm (10 de cada).



USB 2.0



ADJUSTABLE
~20-30x



VGA RESOLUTION
640 x 480



4 WHITE LEDs
NON SWITCHABLE



[EarScope
Product Brochure](#)

Información adicional

- Dino-Lite medical range: Dino-Lite EarScope
- Light/ LED type: White
- Number of LEDs: 4
- LED on/off switchable: Yes
- Infrared filter: IR cut-filter >650 nm
- Diffuser available: No
- Polarizer: No

- Magnification: 20-30x
- Resolution: VGA (640x480)
- Maximum frame rate: 30 fps
- Interface: USB 2.0
- Housing material: Composite/plastic housing
- Magnification lock: No
- Measurement: No
- Calibration: No
- Microtouch sensor: No
- Package contents: Microscope, Carry pouch, 10x 3 mm specula, 10x 4 mm specula, 10x 5 mm specula, Software CD, User manual
- Regulatory approval: Medical Device Class 1 - According to Directive 93/42/EEC modified 2007/47/EC
- Price range: €0,00 - €300,00

• **Details**

- Dino-Lite medical range: Dino-Lite EarScope
- Article code: MEDL3E

• **Lighting**

- Light/ LED type: White
- Number of LEDs: 4
- LED on/off switchable: Yes
- Infrared filter: IR cut-filter >650 nm
- Diffuser available: No
- Polarizer: No

• **Optics**

- Magnification: 20-30x
- Lens type: Glass with Anti-Reflection coating

• **Sensor**

- Sensor type: CMOS
- Resolution: VGA (640x480)
- Maximum frame rate: 30 fps

• **Compatibility**

- Interface: USB 2.0
- Operating system: Windows XP, Vista, 7, 8 & 10, MacOS 10.4 and up
- Included software: DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
- Supported image formats (Windows): BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
- Supported video formats (Windows): WMV, FLV, SWF
- Supported image formats (MacOS): JPEG, PNG
- Supported video formats (MacOS): MOV

- Imaging standards: DirectShow, UVC

- **Housing**

- Housing material: Composite/plastic housing
- Magnification lock: No
- Dimensions: 10.5 cm (L) x 3.2 cm (D)
- Weight: 95 g
- Cable length: 1.8 m

- **Features**

- Measurement: No
- Calibration: No
- Microtouch sensor: No

- **Information**

- Package contents:
Microscope, Carry pouch, 10x 3 mm specula, 10x 4 mm specula, 10x 5 mm specula, Software CD, User manual
- Warranty information: 2 years European warranty
- Regulatory approval: Medical Device Class 1 - According to Directive 93/42/EEC modified 2007/47/EC
- Price range: €0,00 - €300,00

PDF generated on: 27-09-2018