

MS12C



Apoyo de alcance de visión variable con dos adaptadores de alcance de visión variable, transparente. No es compatible con los modelos de polarización y los modelos de la serie AD Pro 2. El MS12C incluye una unidad base con dos adaptadores incluidos que ofrecen tres "manos libres" con opciones de distancia de trabajo fija y de aumento.

Ideal para coleccionistas, aficionados, artistas gráficos, así como para usos en el hogar o en las aulas. Los microscopios digitales Dino-Lite simplemente se deslizan en la base o en uno de los adaptadores. El diseño de la base en acrílico transparente permite una sencilla accesibilidad de las herramientas y de la luz hacia el objeto observado.

Información adicional

- Accessories range: Basic stands
- Dimensions: 8.2cm (W) x 8.2cm (L) x 7.6cm (H)
- Weight: 54 g
- ESD Safe: No
- Compatible with: All Dino-Lite models, except models in AD series and models with polarizer
- Package contents: Stand + 2 adapters
- Price range: €0 - €50,00

• Details

- Accessory range: Basic stands

- Price range: €0 - €50,00
- ESD-Safe: No
- Dimensions: 8.2cm (W) x 8.2cm (L) x 7.6cm (H)
- Weight: 54 g
- Compatible with: All Dino-Lite models, except models in AD series and models with polarizer
- Package contents: Stand + 2 adapters
- Warranty: 1 year European warranty
- Accessory type: Basic stands
- Price range (excl. VAT): €0 - €50,00
- ESD-safe: No

PDF generated on: 24-10-2018