

RK-02



RK-02 es una solución duradera de soporte con cuello de cisne, diseñada para un uso flexible y versátil con un microscopio digital Dino-Lite. La abrazadera en C del RK-02 se ajusta a prácticamente cualquier superficie de hasta 7 cm (~3 pulgadas) de grosor y se puede montar horizontal o verticalmente. Con este soporte robusto, podrá colocar fácilmente el Dino-Lite doblando el brazo de 60 cm de largo y ajustando el cabezal giratorio. La base del RK-02 también se puede girar, lo que proporciona incluso mayor libertad y se puede fijar con el cierre de desenganche rápido para usar sobre superficies irregulares. El brazo con cuello de cisne RK-02 es la opción ideal para usar de manera flexible cualquiera de los modelos Dino-lite, excepto los modelos de grandes aumentos.

Información adicional

- Accessories range: Professional stands
- Dimensions: Gooseneck 60cm (L) x 1.5cm (D)
- Weight: 1.5 kg
- ESD Safe: No
- Compatible with: All Dino-Lite models
- Package contents: Stand, HD-P1 holder
- Price range: €50,00 - €150,00

• Details

- Accessory range: Professional stands
- Price range: €50,00 - €150,00

- ESD-Safe: No
- Dimensions: Gooseneck 60cm (L) x 1.5cm (D)
- Weight: 1.5 kg
- Compatible with: All Dino-Lite models
- Package contents: Stand, HD-P1 holder
- Warranty: 1 year European warranty
- Accessory type: Professional stands
- Price range (excl. VAT): €50,00 - €150,00
- ESD-safe: No

PDF generated on: 24-10-2018