

PodoScope



Gracias al PodoScope, se pueden documentar cambios en la piel del pie y en el propio pie realizando fotografías y haciendo grabaciones también. Estas imágenes se pueden guardar en el expediente del paciente.



USB 2.0



MEASUREMENT
FUNCTIONALITY



ADJUSTABLE
~10-50x



MAGNIFICATION
LOCK



1.3 MEGAPIXEL
1280 x 1024



8 WHITE LEDs
SWITCHABLE



POLARIZER
ANTI-REFLECTION



EXCHANGEABLE
CAP

El podólogo y el paciente pueden ver las imágenes en tiempo real en un ordenador portátil, en una tableta Windows o en un ordenador normal. No es necesaria una fuente de luz externa en la investigación, puesto que el PodoScope cuenta con luces LED incorporadas (ocho luces LED brillantes). La iluminación, así como el propio PodoScope, reciben la energía a través del puerto USB, de forma que el PodoScope siempre está listo y no depende de baterías.



[PodoScope](#)
[Product Brochure](#)

Información adicional

- Dino-Lite medical range: Dino-Lite PodoScope

- Light/ LED type: White
- Number of LEDs: 8
- LED on/off switchable: Yes
- Infrared filter: IR cut-filter >650 nm
- Diffuser available: No
- Polarizer: Yes
- Magnification: 10-50x
- Resolution: 1.3 Megapixel (1280x1024)
- Maximum frame rate: 30 fps
- Interface: USB 2.0
- Housing material: Composite/plastic housing
- Magnification lock: Yes
- Measurement: Yes
- Calibration: Yes
- Microtouch sensor: Yes
- Package contents: Microscope, Carry pouch, Software CD, Calibration target, User manual, N3C-C- Closed cap, N3C-L- Long cap
- Regulatory approval: Medical Device Class 1 - According to Directive 93/42/EEC modified 2007/47/EC
- Price range: €300,00 - €500,00

• **Details**

- Dino-Lite medical range: Dino-Lite PodoScope
- Article code: MEDL4FW

• **Lighting**

- Light/ LED type: White
- Number of LEDs: 8
- LED on/off switchable: Yes
- Infrared filter: IR cut-filter >650 nm
- Diffuser available: No
- Polarizer: Yes, linear

• **Optics**

- Magnification: 10-50x
- Lens type: Glass with Anti-Reflection coating

• **Sensor**

- Sensor type: CMOS
- Resolution: 1.3 Megapixel (1280x1024)
- Maximum frame rate: 30 fps

• **Compatibility**

- Interface: USB 2.0
- Operating system: Windows XP, Vista, 7, 8 & 10, MacOS 10.4 and up

- Included software: DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
- Supported image formats (Windows): BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
- Supported video formats (Windows): WMV, FLV, SWF
- Supported image formats (MacOS): JPEG, PNG
- Supported video formats (MacOS): MOV
- Imaging standards: DirectShow, UVC

- **Housing**

- Housing material: Composite/plastic housing
- Magnification lock: Yes
- Dimensions: 10.5 cm (L) x 3.2 cm (D)
- Weight: 105 g
- Cable length: 1.8 m

- **Features**

- Measurement: Yes
- Calibration: Yes
- Microtouch sensor: Yes

- **Information**

- Package contents:
Microscope, Carry pouch, Software CD, Calibration target, User manual, N3C-C- Closed cap, N3C-L- Long cap
- Warranty information: 2 years European warranty
- Regulatory approval: Medical Device Class 1 - According to Directive 93/42/EEC modified 2007/47/EC
- Price range: €300,00 - €500,00

PDF generated on: 27-09-2018